

Anti-Loxl2 Monoclonal Antibody For Wb, Elisa - Q9Y4K0

Item Number: CM0004929



Introduction

Rabbit monoclonal antibody against human LOXL2 for WB and ELISA applications. Targets lysyl oxidase homolog 2, involved in chromatin remodeling, EMT, and extracellular matrix cross-linking. Cross-reacts with human. ~87kDa.

[Learn More](#)

Parameter	Value
Product Name	LOXL2 Rabbit mAb
Remarks/Alias	LOR; LOR2; WS9-14; LOXL2
Species	Human
Gene ID (Human)	4017
Gene ID	4017
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 675-774 of human LOXL2 (Q9Y4K0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y4K0
Protein Weight	87kDa
Shipping	Ice bag