

Anti-Loxl2 Polyclonal Antibody For Wb, Elisa - Q9Y4K0

Item Number: CM0016999



Introduction

Rabbit polyclonal antibody targeting human LOXL2, validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying lysyl oxidase homolog 2 in EMT, chromatin regulation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	LOXL2 Rabbit pAb
Remarks/Alias	LOR; LOR2; WS9-14; LOXL2
Species	Human
Gene ID	4017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-290 of human LOXL2 (NP_002309.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4K0
Protein Weight	87kDa
Shipping	Ice bag