


Anti-Lpar2 (Edg4, Lpa2) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hbw0

Item Number: CM0018581



Introduction

Rabbit polyclonal antibody against human LPAR2 (EDG4, LPA2) for WB and ELISA. Target is a lysophosphatidic acid receptor involved in G protein-coupled signaling and diverse cellular activities.

[Learn More](#)

Parameter	Value
Product Name	LPAR2 Rabbit pAb
Remarks/Alias	EDG4; LPA2; EDG-4; LPA-2; LPAR2
Species	Human
GeneID (human)	9170
GeneID	9170
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 250-350 of human LPAR2 (NP_004711.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HBW0
Protein Weight	39kDa
Shipping	Ice bag