

Anti-Lpar4 Polyclonal Antibody For Wb, Elisa - Q99677

Item Number: CM0021292

Introduction

Rabbit polyclonal antibody targeting human LPAR4/LPA4 receptor (Q99677). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for GPCR and lysophosphatidic acid signaling research.

[Learn More](#)

Parameter	Value
Product Name	LPAR4 Rabbit pAb
Remarks/Alias	LPA4; P2Y9; GPR23; P2RY9; P2Y5-LIKE; LPAR4
Species	Human
GeneID (Human)	2846
GeneID	2846
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 311-370 of human LPAR4 (Q99677).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99677
Protein Weight	42kDa
Shipping	Ice bag