

Anti-Lpcat1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nf37

Item Number: CM0030249

Introduction

Rabbit polyclonal antibody against human LPCAT1, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for phospholipid metabolism and platelet-activating factor biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	LPCAT1 Rabbit pAb
Remarks/Alias	AYTL2; lpcat; AGPAT9; LPLAT8; PFAAP3; AGPAT10; LPCAT-1; lysoPAFAT; LPCAT1
Species	Human
GeneID	79888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 325-534 of human LPCAT1 (NP_079106.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NF37
Protein Weight	59kDa
Shipping	Ice bag