

Anti-Lpp Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q93052

Item Number: CM0002072

Introduction

Anti-LPP rabbit monoclonal antibody for WB, IHC-P, IF/ICC, ELISA. Recognizes human and mouse lipoma-preferred partner, a scaffold protein in cell adhesion and signaling.

[Learn More](#)

Parameter	Value
Product Name	LPP Rabbit mAb
Remarks/Alias	LPP; lipoma-preferred partner
Species	Human
Gene ID (Human)	4026
Gene ID	4026
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LPP (Q93052).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q93052
Protein Weight	66kDa
Shipping	Ice bag