

## Anti-Lrat Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95237

Item Number: CM0012900

### Introduction

Anti-LRAT rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat LRAT. Suitable for retinoid metabolism and vision research. CamelBio offers this as part of its comprehensive IVD raw material platform.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | LRAT Rabbit pAb                                                                                                  |
| Remarks/Alias    | LCA14; LRAT                                                                                                      |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 9227                                                                                                             |
| GeneID           | 9227                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human LRAT (NP_004735.2). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | O95237                                                                                                           |
| Protein Weight   | 26kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |