


Anti-Lrg1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P02750

Item Number: CM0025741



Introduction

Anti-LRG1 Rabbit Polyclonal Antibody validated for Western blot, immunofluorescence, and ELISA. Detects leucine-rich alpha-2-glycoprotein in human, mouse, rat samples. Suitable for plasma protein studies and biomarker discovery.

[Learn More](#)

Parameter	Value
Product Name	LRG1 Rabbit pAb
Remarks/Alias	LRG; HMFT1766; LRG1
Species	Human
Gene ID (Human)	116844
Gene ID	116844
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-347 of human LRG1 (NP_443204.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02750
Protein Weight	38kDa
Shipping	Ice bag