


Anti-Lrig1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Ja1

Item Number: CM0018091



Introduction

Rabbit polyclonal antibody against human LRIG1, validated for Western blot and ELISA. Detects endogenous LRIG1 in human and mouse samples. Negative regulator of receptor tyrosine kinases.

[Learn More](#)

Parameter	Value
Product Name	LRIG1 Rabbit pAb
Remarks/Alias	LIG1; LIG-1; LRIG1
Species	Human
Gene ID	26018
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 816-1093 of human LRIG1 (NP_056356.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96JA1
Protein Weight	119kDa
Shipping	Ice bag