


Anti-Lrp1 Rabbit Monoclonal Antibody For Wb, Elisa - Q07954

Item Number: CM0003893



Introduction

Rabbit monoclonal antibody against human LRP1 (505 kDa), validated for WB and ELISA. Detects human, mouse, and rat LRP1. Ideal for research on lipid metabolism, endocytosis, and neurodegenerative diseases.

[Learn More](#)

Parameter	Value
Product Name	LRP1 Rabbit mAb
Remarks/Alias	APR; KPA; LRP; A2MR; CD91; APOER; LRP1A; TGFBRS; IGFBP3R; IGFBP-3R; IGFBP3R1; LRP1
Species	Human
GeneID	4035
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 4400-4500 of human LRP1 (Q07954).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07954
Protein Weight	505kDa
Shipping	Ice bag