

Anti-Lrp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q07954

Item Number: CM0013986



Introduction

CamelBio's LRP1 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat LRP1, a 505 kDa endocytic receptor involved in lipid homeostasis, tau clearance, and signaling.

[Learn More](#)

Parameter	Value
Product Name	LRP1 Rabbit pAb
Remarks/Alias	APR; KPA; LRP; A2MR; CD91; APOER; LRP1A; TGFBR5; IGFBP3R; IGFBP-3R; IGFBP3R1; LRP1
Species	Human
GeneID (Human)	4035
GeneID	4035
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 4401-4500 of LRP1 (NP_002323.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07954
Protein Weight	505kDa
Shipping	Ice bag