


Anti-Lrp2 (Megalin) Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P - P98164

Item Number: CM0028612



Introduction

Rabbit polyclonal antibody against human LRP2 (megalin/gp330). Validated for Western blot, IHC-P, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for endocytosis and renal/neural research.

[Learn More](#)

Parameter	Value
Product Name	LRP2 Rabbit pAb
Remarks/Alias	DBS; GP330; LRP-2; LRP2
Species	Human
GeneID (Human)	4036
GeneID	4036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 4446-4655 of human LRP2 (NP_004516.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98164
Protein Weight	522kDa
Shipping	Ice bag