

Anti-Lrp4 Rabbit Polyclonal Antibody For If-P, Elisa - 075096

Item Number: CM0019802



Introduction

LRP4 Rabbit Polyclonal Antibody validated for IF-P and ELISA in human samples. Targets LRP4, a receptor involved in bone formation and neuromuscular junction. UniProt: O75096.

[Learn More](#)

Parameter	Value
Product Name	LRP4 Rabbit pAb
Remarks/Alias	CLSS; CMS17; LRP-4; LRP10; MEGF7; SOST2; LRP4
Species	Human
Gene ID (Human)	4038
Gene ID	4038
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1621-1720 of human LRP4 (NP_002325.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human
SWISS	O75096
Protein Weight	212kDa
Shipping	Ice bag