

Anti-Lrp5 Rabbit Monoclonal Antibody For Wb, Elisa - 075197

Item Number: CM0003545

Introduction

LRP5 Rabbit mAb (O75197) for WB and ELISA. Validated in human and mouse, targets the Wnt coreceptor LRP5 involved in bone development and carcinogenesis. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	LRP5 Rabbit mAb
Remarks/Alias	HBM; LR3; OPS; EVR1; EVR4; LRP7; OPPG; BMND1; LRP-5; LRP-7; OPTA1; PCLD4; VBCH2; LRP5
Species	Human
Gene ID (Human)	4041
Gene ID	4041
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1400-1500 of human LRP5 (NP_002326.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75197
Protein Weight	179kDa
Shipping	Ice bag