

Anti-Lrp5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75197

Item Number: CM0012941

Introduction

Rabbit polyclonal antibody against human LRP5, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Coreceptor in Wnt signaling, bone development, and retinal vasculature.

[Learn More](#)

Parameter	Value
Product Name	LRP5 Rabbit pAb
Remarks/Alias	HBM; LR3; OPS; EVR1; EVR4; LRP7; OPPG; BMND1; LRP-5; LRP-7; OPTA1; PCLD4; VBCH2; LRP5
Species	Human
GeneID (Human)	4041
GeneID	4041
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75197
Protein Weight	179kDa
Shipping	Ice bag