


Anti-Lrp6 Rabbit Polyclonal Antibody For If/Icc, Elisa - 075581

Item Number: CM0017653



Introduction

Rabbit polyclonal antibody against human LRP6, a Wnt coreceptor involved in beta-catenin signaling. Suitable for IF/ICC and ELISA applications. Validated for human samples. UniProt 075581.

[Learn More](#)

Parameter	Value
Product Name	LRP6 Rabbit pAb
Remarks/Alias	ADCAD2; STHAG7; LRP6
Species	Human
GeneID (Human)	4040
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-150 of human LRP6 (NP_002327.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	075581
Protein Weight	180kDa
Shipping	Ice bag