


Anti-Lrp6 Rabbit Polyclonal Antibody For Wb, Elisa - 075581

Item Number: CM0029355



Introduction

This rabbit polyclonal antibody targets human LRP6, a Wnt signaling coreceptor, validated for Western blot (WB) and ELISA, with cross-reactivity to mouse and rat. Ideal for studying beta-catenin signaling and related processes.

[Learn More](#)

Parameter	Value
Product Name	LRP6 Rabbit pAb
Remarks/Alias	ADCAD2; STHAG7; LRP6
Species	Human
Gene ID	4040
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human LRP6 (NP_002327.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	075581
Protein Weight	180kDa
Shipping	Ice bag