

## Anti-Lrpap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30533

Item Number: CM0017290



### Introduction

Rabbit polyclonal antibody targeting human LRPAP1 (RAP, alpha-2-MRAP). Suitable for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse, rat. Functions as molecular chaperone for LDL receptor-related proteins.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | LRPAP1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | RAP; MRAP; A2RAP; HBP44; MYP23; A2MRAP; alpha-2-MRAP; LRPAP1                                                        |
| Species          | Human                                                                                                               |
| GeneID (human)   | 4043                                                                                                                |
| GeneID           | 4043                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 35-357 of human LRPAP1 (NP_002328.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P30533                                                                                                              |
| Protein Weight   | 41kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |