


Anti-Lrpprc Rabbit Polyclonal Antibody (Wb/Ihc-P/If-Icc/Elisa) - P42704

Item Number: CM0026474



Introduction

Rabbit polyclonal antibody against LRPPRC, intended for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for mitochondrial RNA metabolism and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	LRPPRC Rabbit pAb
Remarks/Alias	LSFC; GP130; LRP130; MC4DN5; CLONE-23970; LRPPRC
Species	Human
Gene ID (Human)	10128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1041-1394 of human LRPPRC (NP_573566.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42704
Protein Weight	158kDa
Shipping	Ice bag