

Anti-Lrrc14 Rabbit Polyclonal Antibody For Wb, Elisa - Q15048

Item Number: CM0022522

Introduction

Anti-LRRC14 Rabbit Polyclonal Antibody designed for WB and ELISA applications. LRRC14 negatively regulates Toll-like receptor-mediated NF-kappa-B signaling. Reacts with human LRRC14. Suitable for innate immunity research. UniProt Q15048.

[Learn More](#)

Parameter	Value
Product Name	LRRC14 Rabbit pAb
Remarks/Alias	LRRC14A; LRRC14
Species	Human
Gene ID	9684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human LRRC14 (NP_055480.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15048
Protein Weight	55kDa
Shipping	Ice bag