

Anti-Lrrc4 Polyclonal Antibody For Wb, Elisa - Q9Hbw1

Item Number: CM0013487



Introduction

Rabbit polyclonal antibody against human LRRC4 (NGL-2), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience research and glioma studies.

[Learn More](#)

Parameter	Value
Product Name	LRRC4 Rabbit pAb
Remarks/Alias	NAG14; NGL-2; LRRC4
Species	Human
Gene ID (Human)	64101
Gene ID	64101
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 544-653 of human LRRC4 (NP_071426.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HBW1
Protein Weight	73kDa
Shipping	Ice bag