

Anti-Lrrc61 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bv99

Item Number: CM0017623

Introduction

Anti-LRRC61 rabbit polyclonal antibody raised against full-length recombinant human LRRC61. Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for studying leucine-rich repeat-containing proteins.

[Learn More](#)

Parameter	Specification
Product Name	LRRC61 Rabbit pAb
Remarks/Alias	HSPC295; LRRC61
Species	Human
GeneID (human)	65999
GeneID	65999
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-259 of human LRRC61 (NP_076431.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BV99
Protein Weight	28kDa
Shipping	Ice bag