


Anti-Lrrc7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Nw7

Item Number: CM0016078



Introduction

Rabbit polyclonal antibody targeting human LRRC7 (Densin-180), suitable for WB, IF/ICC, and ELISA. Validated for reactivity with human, mouse, and rat. Ideal for postsynaptic research applications.

[Learn More](#)

Parameter	Value
Product Name	LRRC7 Rabbit pAb
Remarks/Alias	DENSIN; LRRC7
Species	Human
GeneID (Human)	57554
GeneID	57554
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 730-890 of human LRRC7 (NP_065845.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NW7
Protein Weight	173kDa
Shipping	Ice bag