

Anti-Lrrfip1 Polyclonal Antibody For Wb, Elisa - Q32Mz4

Item Number: CM0019327



Introduction

LRRFIP1 Rabbit Polyclonal Antibody validated for WB and ELISA, reacting with human, mouse, and rat. A valuable tool for studying transcriptional repression, smooth muscle proliferation, and TLR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	LRRFIP1 Rabbit pAb
Remarks/Alias	GCF2; TRIP; FLAP1; GCF-2; FLAP-1; HUF1-1; FLIAP1; LRRFIP1
Species	Human
GeneID	9208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 589-809 of human LRRFIP1 (NP_001131024.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q32MZ4
Protein Weight	89kDa
Shipping	Ice bag