

Anti-Lrrfip2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y608

Item Number: CM0019087

Introduction

LRRFIP2 Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for Wnt signaling and TLR signaling research.

[Learn More](#)

Parameter	Value
Product Name	LRRFIP2 Rabbit pAb
Remarks/Alias	HUFI-2; LRRFIP2
Species	Human
GeneID (Human)	9209
GeneID	9209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-340 of human LRRFIP2 (NP_060194.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y608
Protein Weight	82kDa
Shipping	Ice bag