


Anti-Lrrk1 Polyclonal Antibody For Wb, Elisa - Q38Sd2

Item Number: CM0024055



Introduction

Rabbit polyclonal antibody against human LRRK1 for WB and ELISA. LRRK1 phosphorylates RAB7A, regulating intracellular trafficking and osteoclast maturation, with implications in bone metabolism research. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	LRRK1 Rabbit pAb
Remarks/Alias	OSMD; RIPK6; Roco1; LRRK1
Species	Human
GeneID (Human)	79705
GeneID	79705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human LRRK1 (NP_078928.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q38SD2
Protein Weight	225kDa
Shipping	Ice bag