


Anti-Lrrk1 Rabbit Polyclonal Antibody For Wb And Elisa - Q38Sd2

Item Number: CM0023157



Introduction

Rabbit polyclonal antibody targeting human LRRK1 (Leucine-rich repeat serine/threonine-protein kinase 1), validated for Western blot and ELISA. Suitable for studies on intracellular trafficking and osteoclast regulation. Predicted molecular weight 225 kDa.

[Learn More](#)

Parameter	Value
Product Name	LRRK1 Rabbit pAb
Remarks/Alias	OSMD, RIPK6, Roco1, LRRK1
Species	Human
Gene ID (Human)	79705
Gene ID	79705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human LRRK1 (NP_078928.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q38SD2
Protein Weight	225kDa
Shipping	Ice bag