

Anti-Lrrk2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5S007

Item Number: CM0013080



Introduction

This LRRK2 polyclonal antibody (Cat. No. A0859) is produced in rabbit against a recombinant protein containing amino acids 750-910 of human LRRK2. Validated for WB, IHC-P, and ELISA, it recognizes human, mouse, and rat LRRK2 (~286 kDa). Suitable for autophagy, neurodegeneration, and Parkinson's disease research.

[Learn More](#)

Parameter	Value
Product Name	LRRK2 Rabbit pAb
Remarks/Alias	PARK8; RIPK7; ROCO2; AURA17; DARDARIN; LRRK2
Species	Human
Gene ID	120892
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 750-910 of human LRRK2 (NP_940980.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5S007
Protein Weight	286kDa
Shipping	Ice bag