

Anti-Lrrk2 Rabbit Monoclonal Antibody For Wb, Elisa - Q5S007

Item Number: CM0006895



Introduction

High-quality Anti-LRRK2 Rabbit Monoclonal Antibody validated for WB and ELISA. Targets human LRRK2 (PARK8, Dardarin), a key kinase in Parkinson's disease, autophagy, and neuronal function. Ideal for IVD and research.

[Learn More](#)

Parameter	Value
Product Name	LRRK2 Rabbit mAb
Remarks/Alias	PARK8; RPK7; ROCO2; AURA17; DARDARIN; LRRK2
Species	Human
GeneID (Human)	120892
GeneID	120892
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 100-340 of human LRRK2 (NP_940980.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5S007
Protein Weight	286kDa
Shipping	Ice bag