

Anti-Lrrk2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q5S007

Item Number: CM0022675



Introduction

Anti-LRRK2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA; reacts with Human, Mouse, Rat; targets Leucine-rich repeat serine/threonine-protein kinase 2 (LRRK2), a key regulator in neuronal plasticity and autophagy. Ideal for Parkinson's disease research.

[Learn More](#)

Parameter	Value
Product Name	LRRK2 Rabbit pAb
Remarks/Alias	PARK8; RIPK7; ROCO2; AURA17; DARDARIN; LRRK2
Species	Human
GeneID (Human)	120892
GeneID	120892
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 750-910 of human LRRK2 (NP_940980.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5S007
Protein Weight	286kDa
Shipping	Ice bag