

Anti-Lrsam1 Polyclonal Antibody For Wb And Elisa - Q6Uwe0

Item Number: CM0025361

Introduction

LRSAM1 Rabbit Polyclonal Antibody, validated for Western blot and ELISA. Detects LRSAM1 (E3 ubiquitin ligase involved in autophagy and neurodegeneration) in human, mouse, and rat samples. Ideal for Charcot-Marie-Tooth disease research.

[Learn More](#)

Parameter	Value
Product Name	LRSAM1 Rabbit pAb
Remarks/Alias	TAL; CMT2P; RIFLE; LRSAM1
Species	Human
GeneID (Human)	90678
GeneID	90678
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UWE0
Protein Weight	84kDa
Shipping	Ice bag