

Anti-Lrwd1 Polyclonal Antibody For Wb, Elisa - Q9Ufc0

Item Number: CM0020521

Introduction

LRWD1 rabbit polyclonal antibody validated for WB and ELISA in human, mouse, rat samples. Targets leucine-rich repeat and WD repeat-containing protein 1 (CENP-33/ORCA), involved in DNA replication and chromatin organization.

[Learn More](#)

Parameter	Value
Product Name	Lrwd1 Rabbit pAb
Remarks/Alias	ORCA; CENP-33; Lrwd1
Species	Human
Gene ID (Human)	222229
Gene ID	222229
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human Lrwd1 (NP_690852.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UFC0
Protein Weight	71kDa
Shipping	Ice bag