


Anti-Lsd2/Kdm1B/Aof1 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Nb78

Item Number: CM0003641



Introduction

Rabbit monoclonal antibody against human LSD2/KDM1B/AOF1, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for epigenetics research studying histone demethylation and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	LSD2/KDM1B/AOF1 Rabbit mAb
Remarks/Alias	AOF1; LSD2; C6orf193; LSD2/KDM1B/AOF1
Species	Human
Gene ID (Human)	221656
Gene ID	221656
Immunogen	A synthetic peptide corresponding to amino acids 1-100 of human LSD2/KDM1B/AOF1 (NP_694587.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NB78
Protein Weight	92 kDa
Shipping	Ice bag