


Anti-Lsh/Hells Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q9NrZ9

Item Number: CM0007580



Introduction

Rabbit monoclonal antibody against human Lsh/HELLS (Q9NRZ9), validated for WB, IF, IHC, and ELISA. Cross-reacts with mouse and rat. Suitable for chromatin remodeling and DNA methylation research.

[Learn More](#)

Parameter	Value
Product Name	Lsh/HELLS Rabbit mAb
Remarks/Alias	LSH; ICF4; PASG; SMARCA6; Nbla10143; Lsh/HELLS
Species	Human
Gene ID (Human)	3070
Gene ID	3070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 614-838 of human Lsh/HELLS (NP_060533.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRZ9
Protein Weight	97kDa
Shipping	Ice bag