

Anti-Lsm1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15116

Item Number: CM0017413

Introduction

Rabbit polyclonal antibody targeting human LSM1 (SWISS O15116). Validated for Western blot, IF/ICC, and ELISA. Cross-reacts with human and mouse. Useful for studying mRNA degradation and histone mRNA processing.

[Learn More](#)

Parameter	Value
Product Name	LSM1 Rabbit pAb
Remarks/Alias	CASM; YJL124C; LSM1
Species	Human
Gene ID (Human)	27257
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-133 of human LSM1 (NP_055277.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15116
Protein Weight	15kDa
Shipping	Ice bag