

Anti-Lsm11 Polyclonal Antibody For Wb, Elisa - P83369

Item Number: CM0030096



Introduction

Rabbit polyclonal antibody against human LSM11, a key component of the U7 snRNP complex involved in histone mRNA processing. Validated for WB and ELISA, cross-reacts with human and rat. Ideal for cell cycle and RNA processing studies.

[Learn More](#)

Parameter	Value
Product Name	LSM11 Rabbit pAb
Remarks/Alias	AGS8; LSM11
Species	Human
Gene ID (Human)	134353
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 131-360 of human LSM11 (NP_775762.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P83369
Protein Weight	40kDa
Shipping	Ice bag