

Anti-Lsm12 Rabbit Polyclonal Antibody For Wb, Elisa - Q3Mhd2

Item Number: CM0022100

Introduction

Rabbit polyclonal antibody against human LSM12, validated for WB and ELISA. Cross-reacts with mouse. Suitable for NAADP-mediated Ca²⁺ signaling research. Protein weight 22kDa.

[Learn More](#)

Parameter	Specification
Product Name	LSM12 Rabbit pAb
Remarks/Alias	PNAS-135; LSM12
Species	Human
GeneID (Human)	124801
GeneID	124801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human LSM12 (NP_689557.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q3MHD2
Protein Weight	22kDa
Shipping	Ice bag