

## Anti-Lsm14A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nd56

Item Number: CM0019584



### Introduction

Rabbit polyclonal antibody targeting human LSM14A (RAP55), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Key tool for studying P-bodies, translational repression, and mitotic spindle regulation.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | LSM14A Rabbit pAb                                                                                       |
| Remarks/Alias    | RAP55; FAM61A; RAP55A; C19orf13; LSM14A                                                                 |
| Species          | Human                                                                                                   |
| GeneID (Human)   | 26065                                                                                                   |
| GeneID           | 26065                                                                                                   |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human LSM14A (NP_001107565.1) |
| Source           | Rabbit                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                       |
| SWISS            | Q8ND56                                                                                                  |
| Protein Weight   | 51kDa                                                                                                   |
| Shipping         | Ice bag                                                                                                 |