


Anti-Lsm2 Polyclonal Antibody For Wb, Elisa - Q9Y333

Item Number: CM0027575



Introduction

Rabbit polyclonal antibody targeting human LSM2, validated for WB and ELISA. Detects an 11 kDa protein, with reactivity confirmed in human samples. Optimized for pre-mRNA splicing and spliceosome research applications.

[Learn More](#)

Parameter	Value
Product Name	LSM2 Rabbit pAb
Remarks/Alias	G7B; snRNP; C6orf28; YBL026W; LSM2
Species	Human
GeneID (Human)	57819
GeneID	57819
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y333
Protein Weight	11kDa
Shipping	Ice bag