

Anti-Lsm4 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Y4Z0

Item Number: CM0028461



Introduction

Rabbit polyclonal antibody against LSM4, a U6 snRNA-associated Sm-like protein involved in pre-mRNA splicing. Validated for Western blot (WB), immunoprecipitation (IP), and ELISA; cross-reacts with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	LSM4 Rabbit pAb
Remarks/Alias	GRP; YER112W; LSM4
Species	Human
GeneID (Human)	25804
GeneID	25804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-139 of human LSM4 (NP_036453.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4Z0
Protein Weight	15kDa
Shipping	Ice bag