

Anti-Lsp1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P33241

Item Number: CM0005927



Introduction

LSP1 Rabbit mAb (P33241) validated for WB, IF/ICC, IF-P, ELISA in human, mouse, rat. Specific to lymphocyte-specific protein 1, ideal for neutrophil activation and chemotaxis research.

[Learn More](#)

Parameter	Value
Product Name	LSP1 Rabbit mAb
Remarks/Alias	WP34; pp52; LSP1
Species	Human
GeneID	4046
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33241
Protein Weight	37kDa
Shipping	Ice bag