

Anti-Lsp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P33241

Item Number: CM0029547



Introduction

Rabbit polyclonal antibody against human LSP1, validated for WB, IHC-P, and ELISA. Cross-reacts with human and mouse. Useful for studying neutrophil activation and chemotaxis.

[Learn More](#)

Parameter	Value
Product Name	LSP1 Rabbit pAb
Remarks/Alias	WP34; pp52; LSP1
Species	Human
GeneID (human)	4046
GeneID	4046
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human LSP1 (NP_002330.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P33241
Protein Weight	37kDa
Shipping	Ice bag