

Anti-Lsr Polyclonal Antibody For Wb And Elisa - Q86X29

Item Number: CM0029455

Introduction

Rabbit polyclonal antibody targeting human LSR (lipolysis-stimulated lipoprotein receptor), validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studies of lipoprotein clearance and tight junction biology.

[Learn More](#)

Parameter	Value
Product Name	LSR Rabbit pAb
Remarks/Alias	ILDR3; LISCH7; LSR
Species	Human
GeneID (Human)	51599
GeneID	51599
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-240 of human LSR (NP_057009.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86X29
Protein Weight	71kDa
Shipping	Ice bag