


Anti-Lta Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01374

Item Number: CM0016616



Introduction

High-quality rabbit polyclonal antibody against human LTA (Lymphotoxin-alpha), validated for WB, IHC-P, and ELISA. Reacts with human and mouse. Ideal for cytokine and tumor necrosis research.

[Learn More](#)

Parameter	Value
Product Name	LTA Rabbit pAb
Remarks/Alias	LT; TNFB; TNFSF1; TNLG1E; LTA
Species	Human
Gene ID (Human)	4049
Gene ID	4049
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-205 of human LTA (NP_000586.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01374
Protein Weight	22kDa
Shipping	Ice bag