


Anti-Lta4H Monoclonal Antibody For Wb - P09960

Item Number: CM0009961



Introduction

Rabbit monoclonal antibody against human LTA4H (leukotriene A4 hydrolase) for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for inflammation and lipid mediator research.

[Learn More](#)

Parameter	Specification
Product Name	LTA4H Rabbit mAb
Remarks/Alias	LTA4H; leukotriene A4 hydrolase
Species	Human
Gene ID (Human)	4048
Gene ID	4048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 501-600 of human LTA4H (P09960)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09960
Protein Weight	69kDa
Shipping	Ice bag