


Anti-Lta4H Polyclonal Antibody For Wb And Elisa - P09960

Item Number: CM0020336



Introduction

Rabbit polyclonal antibody against human LTA4H, validated for WB and ELISA. Cross-reacts with mouse. LTA4H is a bifunctional enzyme involved in leukotriene B4 biosynthesis and inflammation resolution. Ideal for inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	LTA4H Rabbit pAb
Remarks/Alias	LTA4H
Species	Human
GeneID (Human)	4048
GeneID	4048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 312-611 of human LTA4H (NP_000886.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09960
Protein Weight	69kDa
Shipping	Ice bag