


Anti-Ltb4R Rabbit Polyclonal Antibody For Wb, Elisa - Q15722

Item Number: CM0013066



Introduction

LTB4R Rabbit pAb is a polyclonal antibody raised in rabbit against a synthetic peptide from human LTB4R. Validated for WB and ELISA, it cross-reacts with mouse and rat, enabling research on inflammation and immune response.

[Learn More](#)

Parameter	Value
Product Name	LTB4R Rabbit pAb
Remarks/Alias	BLT1; BLTR; P2Y7; GPR16; LTBR1; P2RY7; CMKRL1; LTB4R1; LTB4R
Species	Human
Gene ID	1241
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 253-352 of human LTB4R (NP_001137391.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q15722
Protein Weight	38kDa
Shipping	Ice bag