

Anti-Ltbp3 Polyclonal Antibody For Wb, Elisa - Q9Ns15

Item Number: CM0013133

Introduction

Anti-LTBP3 rabbit polyclonal antibody for WB and ELISA, reactive with human and rat. LTBP3 is a key regulator of TGF-beta activation, associated with extracellular matrix. Ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	LTBP3 Rabbit pAb
Remarks/Alias	DASS; LTBP2; LTBP-3; STHAG6; pp6425; GPHYSD3; LTBP3
Species	Human
GeneID (Human)	4054
GeneID	4054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1004-1303 of human LTBP3 (NP_001123616.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9NS15
Protein Weight	139kDa
Shipping	Ice bag