

Anti-Ltbr Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36941

Item Number: CM0029136



Introduction

Rabbit polyclonal antibody targeting human LTBR/TNFRSF3. Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for apoptosis and NF-κB signaling studies.

[Learn More](#)

Parameter	Value
Product Name	LTBR Rabbit pAb
Remarks/Alias	TNFCR; TNFR3; D12S370; TNFR-RP; TNFRSF3; TNFR2-RP; LT-BETA-R; TNF-R-III; LTBR
Species	Human
Gene ID (Human)	4055
Gene ID	4055
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 336-435 of human LTBR (NP_002333.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36941
Protein Weight	47kDa
Shipping	Ice bag