

Anti-Ltf Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P02788

Item Number: CM0017467



Introduction

Human LTF rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat lactotransferrin (~78 kDa), a key iron-binding protein in innate immunity and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	LTF Rabbit pAb
Remarks/Alias	LF; HLF2; GIG12; HEL110; LTF
Species	Human
Gene ID	4057
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-352 of human LTF (NP_002334.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02788
Protein Weight	78kDa
Shipping	Ice bag